

# Platelet-Rich Plasma (PRP): Standard Operating Procedure (SOP)

## Document Control

Document Number	PRP-SOP-001
Version	1.1
Effective Date	2026 APRIL
Review Date	2027 APRIL

## Approval

Prepared by: Dr Albert	Reviewed by: Dr Amy	Approved by:
Name/Signature: Albert	Name/Signature: Amy	Name/Signature:
Date: 2026 APRIL	Date: 2026 APRIL	Date:

## 1. Purpose

To ensure proper preparation, storage, handling, and appropriate laboratory use of Platelet-Rich Plasma (PRP) for transfusion.

---

## 2. Scope

This SOP applies to all laboratory personnel involved in:

- Preparation of PRP from whole blood
  - Storage and monitoring
  - Issue and documentation
- 

## 3. Responsibility

- **Pathologist / Transfusion Specialist:** Overall supervision and approval
  - **Medical Laboratory Technologist:** Preparation, validation, and release
  - **Laboratory Technician:** Processing, storage, labeling, and documentation
- 

## 4. Definitions

- **Platelet-Rich Plasma (PRP):** A whole-blood-derived platelet component prepared by the PRP method (light centrifugation of whole blood), used as the available platelet product in this laboratory
- 

## 5. Procedure

### 5.1 Preparation of PRP

- PRP is prepared from fresh whole blood **within 6–8 hours of collection**
- Use **light (soft) centrifugation** [light spin using  $2000 \times g$  for 3 minutes (plus deceleration time) with a temperature setting of  $20^{\circ} C$ ] to separate PRP from red cells
  - To calculate RPM for the centrifuge, the following formula is used:
  - **$RCF = 28.38 \times R \times (RPM/1000)^2$  or  $RPM = \sqrt{RCF / (28.38 \times R)} \times 1000$**

- [RCF=relative centrifugal force ( $\times$  g); R=radius in inches; RPM=revolutions per minute]
  - Transfer the supernatant plasma (PRP) into a satellite bag using sterile technique
  - Avoid contamination and excessive agitation during preparation
- 

## 5.2 Storage of PRP

- Store PRP at **20–24°C**
  - Maintain continuous **gentle agitation** to preserve platelet function
  - **Do NOT refrigerate** (cold temperature reduces platelet viability)
  - Shelf life:
    - Up to **72 hours** for PRP prepared by the whole-blood PRP method, unless a validated and approved system allows extended storage
  - Monitor storage temperature regularly
  - Units must be properly labeled (donor ID, date/time, expiry)
- 

## 5.3 Inspection Before Issue

Each PRP unit must be checked for:

- Swirling (indicator of viable platelets)
- Absence of clots
- No discoloration or contamination
- Intact bag and labeling

Units failing inspection must be discarded.

---

## 5.4 Clinical Indications

PRP transfusion is indicated for:

- Thrombocytopenia with bleeding
  - Platelet count  $<10 \times 10^9/L$  (prophylactic transfusion, except in ITP/TTP/HUS)
  - Platelet count  $<50 \times 10^9/L$  with active bleeding or invasive procedures
  - Platelet count  $<70 \times 10^9/L$  for surgery with anticipated major blood loss
  - Platelet count  $<100 \times 10^9/L$  for critical-site surgery (e.g., brain, eye)
  - Platelet dysfunction with bleeding
- 

## 5.5 Not Indicated For

- Stable patients without bleeding and adequate platelet count
  - As a volume expander
  - Thrombocytopenia due to platelet destruction (unless clinically justified)
- 

## 5.6 Issue and Handling

Before issue:

- Verify:
  - Patient identity (**two identifiers**)
  - ABO compatibility is preferable but not mandatory; **crossmatching is not routinely required for PRP/platelet products**

- Record:
  - Unit number
  - Date/time of issue
  - Patient details

After issue:

- Transfuse as soon as possible after issue and **complete transfusion within 30–60 minutes** (not exceeding **6 hours from issue** as per local policy)

---

### 5.7 Post-Issue Handling

- PRP should not normally be returned to inventory once issued.
- A PRP unit may only be returned to stock if ALL of the following conditions are met:
  - The unit has not been transfused, entered, or compromised
  - The unit has been maintained at 20–24°C with continuous gentle agitation
  - The unit has been handled within validated institutional procedures
  - The unit is visually acceptable (no clots, contamination, discoloration, or damage)
  - The unit identity and labeling are intact and traceable
  - Return is reviewed and approved by authorized laboratory personnel
- Discard the PRP unit if any of the following occur:
  - Storage temperature or agitation conditions are not maintained within acceptable limits
  - Unit is expired
  - Evidence of contamination, clots, leakage, or abnormal appearance
  - Unit identity or traceability is compromised

---

### 6. Quality Control and Assurance

- Monitor:
  - Storage temperature (20–24°C)
  - Agitation system function
- Maintain:
  - Traceability of each unit
  - Complete documentation
- Perform:
  - Equipment calibration
- Participate in:
  - Quality assurance programs

---

### 7. Safety Precautions

- Use PPE at all times
- Treat all blood products as infectious
- Follow biosafety and waste disposal protocols

---

### 8. Documentation

Maintain records of:

- Donor information
  - Preparation details
  - Storage logs
  - Issue records
  - Discard records
- 

## 9. References

- AABB Standards for Blood Banks and Transfusion Services
  - AABB Technical Manual
  - Guidelines for Blood Transfusion Service, MOHS, Myanmar, 2018
- 

## Audit Statement

This SOP is based on AABB standards and aligned with MOHS guidelines where applicable.